

**skin condition**

Balanced to mature or prematurely-aging skin.

**description**

Pro Microneedling is an advanced treatment that delivers innovative techniques and powerful, professional-grade actives to brighten the complexion, refine the look of pores, reduce the appearance of dark spots, and diminish the appearance of lines and wrinkles.

**benefits**

- Brightens complexion
- Reduces the appearance of dark spots
- Diminishes the appearance of lines and wrinkles
- Refines the look of pores

**featured professional products**

- Pro Power Peel
- Hyaluronic Acid IonActive
- PowerBright IonActive Serum
- Pro Restore



# pro microneedling

**how it works****1. Preparation**

Skin resurfacing with customized Pro Power Peel

**2. Treatment**

Microneedling with professional serums

**3. Recovery**

In-treatment room recovery along with post-treatment care and home care

**treatment protocol****double cleanse**

1. PreCleanse
2. Perform Face Mapping
3. Skin specific cleanser

**resurfacing**

1. Apply One Step Prep.
2. Apply specific Pro Power Peel (1-2 layers).
3. Apply Neutralizing Solution after Pro Power Peel and remove as directed.
  - a. Power Down: Use Exfoliant Accelerator 35 or omit resurfacing.

**deep treatment**

Optional: Apply numbing cream and allow to process per product guidelines. Remove with one cleanse of Special Cleansing Gel and degrease the skin with One Step prep or Post Extraction Solution.

1. Apply skin specific professional serum.
2. Perform the Dermalogica Microneedling Technique.
  - a. Wipe the skin with damp esthetic wipes after microneedling.
  - b. If including neck and chest, add additional 10 minutes per area.

- Dull Skin: BioLumin-C Pro Serum or Pro Restore at 0.2-0.6mm
- Pore Refining: Niacinamide IonActive at 0.4-0.8mm or Pro Restore 1.0-1.2mm
- Hyperpigmentation: PowerBright IonActive Serum or Pro Restore at 0.2-0.8mm
- Lines and Wrinkles: Pro Restore at 0.4-1.2mm

**dermal layering**

1. Apply Pro Restore.