

skin condition

14 pin cartridge: acne scarring (PIH, PIE, Rolling scars, Box Car scars, and Ice Pick scars)

Metal nano + Silicon nano cartridges: NanoInfusion of professional-grade products to treat dehydration, fine lines, wrinkles, and pigmentation.

description

The Pro Pen Microneedling System is a microneedling device and accessories intended to be used as a treatment to improve the appearance of facial acne scars in adults aged 22 years and older.

It is a multipurpose microneedling and nanoinfusion system with interchangeable cartridges that allow professionals to adjust the degree of treatment and support customized client experiences.



professional only



pro pen

application**microneedling and nanoinfusion**

1. Choose the appropriate cartridge
2. Adjust speed and depth
3. Follow the Dermalogica Technique
 - Apply professional serum or Pro Restore zone by zone
 - Move perpendicular (90°) to skin in slow S-pattern or snake-like movement with horizontal and vertical passes on each area

3 cartridge options allow for customization for every client**• 14-Pin**

FDA-cleared to improve the appearance of facial acne scars.
Stainless-steel needles on 7.86mm round plate
Provides more evenly distributed microchannels
Easily identifiable by its blue color

• Nano Metal

NanoInfusion of professional-grade products to treat dehydration, fine lines, wrinkles, + pigmentation.
Fine cones on 5×5mm stainless-steel plate
Larger plate covers broader areas

• Nano Silicon

NanoInfusion of professional-grade products to treat dehydration, fine lines, wrinkles, + pigmentation.
Fine cones on 2.87×2.87mm silicon plate
Designed for delicate areas, such as eyes and lips